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25.04.2002/CB/GP

Operating Manual Mode d'emploi

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Deutsch

High precision Tubing Pump with planetary drive

IP

| | |
|-------------|---------|
| 4 Channels | ISM 940 |
| 8 Channels | ISM 941 |
| 12 Channels | ISM 942 |
| 16 Channels | ISM 943 |
| 24 Channels | ISM 944 |

IP-N

| | |
|-------------|---------|
| 4 Channels | ISM 945 |
| 8 Channels | ISM 946 |
| 12 Channels | ISM 947 |
| 16 Channels | ISM 948 |
| 24 Channels | ISM 949 |

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Contents

| | |
|-----------------------------|----|
| Pumping | 22 |
| Pumping against pressure | 23 |
| When the pump is not in use | 23 |
| Overload protector | 24 |
| | |
| Tubing | 25 |
| Running-in period | 25 |
| Tubing life | 25 |
| | |
| Analog interface | 26 |
| Settings of DIP-Switch | 29 |
| | |
| Accessories | |
| Foot switch | 30 |
| Spare cassettes | 30 |
| | |
| Flow rates charts | 32 |
| | |
| Maintenance | 35 |
| Repairs | 35 |
| Spare-parts | 36 |
| | |
| Technical specifications | 38 |

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Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry.

As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes

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Safety precautions

- The pump must be switched off when cassettes or tubing are inserted or changed.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Pump-heads consist of revolving parts. Therefore, the pump must not be operated before all cassettes are inserted and fully snapped-in.
- Do not manipulate the pump-head or the cassettes before the pump is switched off and disconnected from the mains supply.

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Safety precautions

- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.
- Only new fuses, according to the specifications stated on page 12 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.

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Safety precautions

- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC SA denies any responsibility.

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Warranty

2 years from date of delivery.

Tubes are regarded as consumables and as such they are not covered by the warranty.

In case of any queries, please contact your ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC SA the warranty terms of the specific manufacturer are valid.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs are charged to the customer.

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Warranty terms

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications

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Product

Contents of the package

- IP or IP-N pump
(type as ordered)
incl. Click'n'go cassettes (POM)
- 1 power cord with an IEC 320
connector (female plug) and a
country specific mains plug
- Operating manual

Please check the package
and its contents for transport
damage. If you find any signs of
damage, please contact your local
ISMATEC® representative
immediately.

➔ Complaints can only be accepted
within 8 days from receipt of the
goods.

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
Rear panel


- 1 Analog interface
Input for:
 - speed control
0–5 V or 0–10 V, and
0–20 mA or 4–20mA
 - Run/Stop
 - rotation direction
- 2 Mains switch (on/off)
- 3 Fuse-holder with voltage
selector 115/230 V
- 4 Mains socket

**⚠ For further details on the
analog interface please refer to
page 26–29.**

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- 1 LED-Test »8.8.8.8.«
- 2 Version der Systemsoftware

Starting the pump

- Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage (page 13).
- Prepare the cassettes (page 15)
- Insert the tubing (page 16)
- Insert the cassettes (page 16)
- Connect the pump tubing
- Connect the pump to the mains and switch it on with the power supply switch

Start-up information

After switching on the power supply switch, the following values are displayed:

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
⚠ Before starting-up
Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ Socket/Power cord
Use exclusively the original power cord supplied with the pump. The socket must be earthed. (protective conductor contact)

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
Voltage setting 115V / 230V Changing the fuses

- ➔ Switch the pump off,
pull out the mains plug.
- 1 Pull out the fuse-holder by opening
it with a small screw-driver (size 0).
 - 2 Take out the voltage selector plate.
Turn it and re-insert it into the fuse-
holder so that the required voltage
rating is facing the window of the
fuse-holder.
 - 3 Insert 2 new fuses
230 V_{AC}: 2 x 250 mA (slow-blow)
115 V_{AC}: 2 x 500 mA (slow-blow)


⚠ Use always 2 slow-blow fuses
of the same type complying with
the local mains voltage.
 - 4 Insert the fuse-holder. The voltage
rating is visible in the window.

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Preparing the cassettes

- 1 ➔ Switch the pump off
- 2 Remove the cassette by slightly pressing the release lever on both sides of the pump.
- 3 The cassettes have slits on both sides. Insert a cassette adaptor into the slit on each side of the cassette.

Press the adaptor firmly against the bottom of the slit.


- ➔ The better you insert the adaptor the easier you can insert the cassettes loaded with the tubing.

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Inserting the tubing


- 1 ➔ Switch the pump off
- 2 Insert the 2-stop tubing with one stop into an adaptor. Stretch the tubing slightly so that it can be inserted with the second stop into the other adaptor.
➔ Ensure that the tube is not twisted.
- 3 Reinsert all cassettes into the pump. Ensure that they snap in on both sides.
➔ The direction of the arrow engraved on the cassette indicates the flow direction providing the least pulsation
- 4 Switch the pump on and start it
(➔ see also information on p. 31)

➔ When the pump is idle, we recommend you to release the tubes from pressure by disengaging all cassettes on one side.

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Operating panel

- 1 Digital LED display
- 2 Increase speed
- 3 – Reduce speed
– settings*
- 4 – MIN key (for standby)
– default*
- 5 – MAX key (ideal for filling or
emptying the tubing)
– ok*
- 6 – Rotation direction
– escape*
- 7 RUN/STOP
- 8 Operating mode
– **PUMP (local operation)**
– REMOTE (control via analog
interface (page 26–29))

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Control keys

The control keys are not active under remote control

a RUN/STOP

- starts the pump at the preset speed
(the rotation direction is briefly indicated with a running dash)
- stops the pump

b Rotation direction

- changes the rotation direction
➡ minus sign = counter-clockwise direction

c ▲ ▼ Keys

- ▲ = increase speed
- ▼ = reduce speed

d MAX

- max. rotation speed
(keep key pressed whilst the pump is running)

e MIN

- The pump runs at 4% of the max. speed (e.g. for stand-by over night)

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
Control keys for the basic settings

- a settings**
Allows access to the menu settings
- b ok**
Confirms a data entry
- c escape**
 - Cancels an entry
 - Quits the menu settings
- d default**
Sets a value to its default value
- e ▲▼ Keys**
 - Select the basic settings
 - Select-value

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Selecting the basic settings

- 1 Push the settings key whilst switching the pump on
- 2 Select adjustable operating settings with the ▲▼ keys
 - FS (foot switch)
 - Speed (rotation speed)
 - Init (default values)
- 3 ok
Confirm your selection
- 4 Select the value by using the ▲▼ key, or
 - press escape for cancelling wrong input
 - certain values can be reset to the default value by using the default key
- 5 ok
Confirm entered value
- 6 escape
Quit menu basic settings

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Basic settings

Foot switch or hand dispenser

- Toggle (default value)
 - Step on foot switch
= pump starts running
 - Press foot switch again
= pump stops
- Direct: the pump runs as long as the foot switch is pressed

Speed

- % (default value)
Rotation speed in percent
Setting: 1.0 – 100.0 %
Resolution: 0.1 %
- **RPM**
Rotation speed in rpm
Setting: 0.4 – 45.0 rpm (IP)
0.11– 11.25 rpm (IP-N)
Resolution
0.1 rpm (IP) / 0.03 rpm (IP-N)

Initializing (default value)

The ok key resets all basic settings to their default values

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Pumping

- 1 Set the required rotation speed (in rpm or %) with the ▲▼ keys
 - 2 Enter the rotation direction with the ◀▶ key
 - 3 Start the pump with the RUN/STOP key
- ➔ The flow rate can also be adjusted whilst the pump is running
 - ➔ For fast filling or emptying the system, maintain pressure on the MAX key whilst the pump is running
 - ➔ For stand-by, e.g. during the night, push the MIN key. The pump runs at 4% of the maximum speed.

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Pumping against pressure

The IP pump can be used for continuous duty at a differential pressure of max. 1.0 bar (smaller tube sizes and/or lever type cassettes allow higher pressures)

In case of any doubts please contact your ISMATEC® agent.

When the pump is idle

When the pump is idle, we recommend you to release the tubing from pressure. Releasing the cassette on the right side is sufficient. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

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Overload protector

The IP pump features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

OL I = Current
OL S = Speed
OL t = Temperature

Whenever this situation occurs, the pump must be switched off immediately.

Before the pump is re-started, it is most important to check the reason for the overload (e.g. too high differential pressure or pressure lever tightened too much, etc.).

➔ Only when the operating conditions have been changed accordingly may the pump be started again.

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Please note

Please also refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and repeatable flow rates are required, we recommend you to run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation!).

In order to reduce wear and tear of the tubing, we recommend users to spray both the tubing and the pump rollers with silicone oil spray (Order No. SC0179). This may improve the service life of the tubing.

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
Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, start

- In normal operation (pin 2 open) pin 3 operates the foot-switch (see operating settings, page 21)
- In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 4, direction

In the open position the pump turns clockwise; when connected to pin 1 (GND) it turns counter-clockwise.

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
External speed control

(0–5V, 0–10V, 0–20mA, 4–20mA)

Input impedance and input range can be selected via a dip-switch inside the pump (see p. 29)

Pin 7, +30 V_{DC}

About +30 V_{DC} are available (max. current 0.2 A).

Pin 9, speed OUT

0–10 V_{DC}, proportionally to the speed

0–45.0 rpm (IP) or

0–11.0 rpm (IP-N).

Pin 10, +5 V_{DC}

About +5 V_{DC} are available. (max. current 0.1 A)

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
Analog interface

Pin 13, speed intern

Depending on how the pump is operated, pin 13 has different functions:

Analog interface not activated

(Normal operation, i.e. pin 2 is open) Pin 13 serves as **auto-start function**.

If pin 13 is connected to pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be on).

Analog interface activated

(Pin 2 on GND)

- Pin 13 open:
The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND:
The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
Settings of DIP-Switch

➔ Make sure that the pump is disconnected from the mains supply.

⚠ The instrument should only be opened by a qualified technician!
Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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| Pins | | | DIP-Switch 1 | DIP-Switch 2 | DIP-Switch 3 | DIP-Switch 4 | DIP-Switch 5 | DIP-Switch 6 |
|-------------------|--------|------|--------------|--------------|--------------|--------------|--------------|--------------|
| Pin 5 speed IN | 0–5V | 1MΩ | OFF* | OFF* | OFF* | | | |
| | 0–10V | 2kΩ | OFF | ON | OFF | | | |
| | 0–20mA | 250Ω | ON | OFF | OFF | | | |
| | 4–20mA | 250Ω | ON | OFF | ON | | | |

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Accessories

Foot switch

Order No. IS 10039
This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

Spare cassettes in POM

CA Click'n'go (since Nov. 2001 with stainless steel pressure spring)
Order No. IS 3710

CA pressure lever*
Order No. IS 0122

Spare cassettes in PVDF

CA pressure lever*
Order No. IS 3820

Adaptor for type CA IS 0123
(order 2 adaptors per cassette)

* The cassettes with pressure lever are available on request. Depending on the application, this type of cassette may provide better results at elevated differential pressure conditions.

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Accessories

Color-coded 2-stop tubing

For reproducible results we recommend you to use tubing from IFMATEC®. We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

⚠ When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you to prime the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes. For applications with a high differential pressure (>1 bar), we recommend the use of cassettes with pressure lever.

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⚠ Flow rates

In the tables listed subsequently the values indicated are only approximate. They were determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

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Please note

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you to use tubing from ISMATEC®

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Flow rates

⚠ In the tables listed subsequently the values indicated are only approximate. They were determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

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Maintenance

Provided the pump is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective pump with detailed description of the failure to your ISMATEC® representative.

Please use the original ISMATEC® packing or a packet of equal quality.

Please give information on defect, date of purchase, serial-no., and model.

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Spare-parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask from your ISMATEC® representative:

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Motor type DC motor

Speed

0.4 – 45.0 rpm (IP)
0.11 – 11.25 rpm (IP-N)
digitally adjustable in steps
of 0.1% or
0.1 rpm (IP)
0.03 rpm (IP-N)

Pump rollers

8 rollers, stainless steel
actively driven

Differential pressure

max. 1.0 bar (14.5 psi)
(For short term use, higher pressures are possible provided that small tubing i.d.'s and/or cassettes with pressure levers are used)

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Please note

We also recommend you to observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Temperature +41 to 104°F

Rel. humidity max. 80%

– not condensing, at normal laboratory conditions

Dimensions/Weight

| Channels | DxWxH (inches) | Weight |
|----------|----------------|---------|
| 4 | 7.1x5.7x5.1 | 10.1 lb |
| 8 | 8.7x5.7x5.1 | 11.2 lb |
| 12 | 10.2x5.7x5.1 | 12.8 lb |
| 16 | 11.8x5.7x5.1 | 14.3 lb |
| 24 | 15.0x5.7x5.1 | 17.4 lb |

Width with adapters: 7.0 inches

CE-compatibility proved according to:

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
Verlangen Sie den
ISMATEC-Katalog.

Ask for the latest
ISMATEC catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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A Unit of IDEX Corporation
Laboratoriumstechnik GmbH