## BioVersal Sample Loops Product Data Sheet



Ultra-high



Face Sealing Technology



PLS Tubing



Learn More at





The BioVersal™ sample loops are 100% bioinert, making them ideal when injecting biological or "sticky" samples. Loops are available in a variety of volumes to enhance application flexibility. No tools are required — these sample loops are supplied with two finger-tight MarvelX™ fittings which are compatible with 10-32 threaded ports. The BioVersal sample loops utilize our unique patented MarvelX technology to auto adjust to various port depths. The outcome — zero dead volume and better chromatographic results.


19,000 psi


FLAT-FACED SEALING CONNECTIONS DEAD VOLUME




















Bioinert Sample Loops					
Part Number	Description	Pressure Rating	Number of Uses	ID μm	Length
UPFP-6005	5 μL Peek Lined SST Sample Loop	1,310 bar (131 MPa, 19,000 psi)	≥ 100	226	125 mm (4.92")
UPFP-6010	10 μL Peek Lined SST Sample Loop	1,310 bar (131 MPa, 19,000 psi)	≥ 100	226	250 mm (9.84")
UPFP-6020	20 μL Peek Lined SST Sample Loop	1,310 bar (131 MPa, 19,000 psi)	≥ 100	320	250 mm (9.84")
UPFP-6050	$50\mu\text{L}$ Peek Lined SST Sample Loop	1,310 bar (131 MPa, 19,000 psi)	≥ 100	320	625 mm (24.60")
UPFP-6100	100 µL Peek Lined SST Sample Loop	1,310 bar (131 MPa, 19,000 psi)	≥ 100	320	1,250 mm (49.21")

Note: All loops use 1/32" PEEK Lined Stainless Steel and include two (2) MarvelX fittings (nut and ferrule).

## **REGULATORY COMPLIANCE**

As of the date of the publication, BioVersal Sample Loops are compliant with current RoHS regulations.



