


## 외경이 다른 튜빙-튜빙 연결 Selection Guide

### 1) 1/16" OD - 1/8" OD 튜빙 연결 시




**2. XP-235**  
(For 1/16" OD Tubing)

**1. P-703-01 (Union Body)**

**3. XP-335**  
(For 1/8" OD Tubing)

### 2) 1/16" OD - 1/32" OD 튜빙 연결 시




(For 1/16" OD Tubing)

(For 1/32" OD Tubing)

**1. P-881**

### 3) 1/16" OD - 360um OD 튜빙 연결 시




(For 1/16" OD Tubing)

(For 360um OD Tubing)

**1. P-882**

4) 1/16" OD - Capillary (<1/32" OD) 튜빙 연결 시



(For 1/16" OD Tubing)

1. P-779

2. 사용할 튜빙 OD Size에 맞는 슬리브선택 (밑의 표 참고)

**NANOTIGHT FEP TUBING SLEEVES, 1/16" OD**

|        | ID             | For Tubing OD Size | Color   | Qty.  |
|--------|----------------|--------------------|---------|-------|
| F-237x | 125µm (.005")  | 70 - 110µm         | Red     | 10-pk |
| F-238x | 180µm (.007")  | 125 - 165µm        | Yellow  | 10-pk |
| F-239x | 215µm (.0085") | 160 - 200µm        | Natural | 10-pk |
| F-240x | 280µm (.011")  | 225 - 265µm        | Blue    | 10-pk |
| F-241x | 330µm (.013")  | 275 - 315µm        | Orange  | 10-pk |
| F-242x | 395µm (.0155") | 340 - 380µm        | Green   | 10-pk |
| F-243x | 455µm (.018")  | 400 - 440µm        | Black   | 10-pk |
| F-244x | 535µm (.021")  | 480 - 520µm        | Natural | 10-pk |
| F-245x | 610µm (.024")  | 555 - 595µm        | Red     | 10-pk |
| F-246x | 685µm (.027")  | 630 - 670µm        | Yellow  | 10-pk |
| F-247x | 840µm (.033")  | 785 - 825µm        | Green   | 10-pk |
| F-252x | 1.07mm (.042") | 1mm                | Purple  | 10-pk |

5. 1/16" OD - 1mm OD 튜빙 연결 시



(For 1/16" OD Tubing)

1. P-702

2. F-252  
(Sleeve, For 1mm OD Tubing)